

HOW TO CLEAN IN-THE-EAR HEARING AIDS

TIPS

- Always wash your hands well before cleaning your hearing aids.
- Leave hearing aids out while showering, washing your face, and washing your hair.
- Always use hair products **BEFORE** putting hearing aids on.
- Clean devices at the end of each day.
DO NOT USE ALCOHOL OR ALCOHOL WIPES.

DAILY CLEANING

1. Examine the hearing aid for any debris and remove with a soft brush or a dry cloth. Be sure to brush microphone covers to remove any debris. Make sure you are brushing with the device held downward, so any debris will fall out of the device and not into it.
2. Wipe the hearing aid down with a dry cloth.
3. Gently remove any earwax that is blocking the sound hole.

DISPOSABLE BATTERIES

If you have a hearing aid with a disposable battery, remove the battery completely once a week and use a cotton swab to gently clean the area where the battery sits. This will help prevent moisture build-up and corrosion.

Continued on back side >



THE
SPEECH & HEARING
CENTER

Celebrating 70 Years - 1953-2023

423.622.6900

www.SpeechHearing.com

HOW TO CLEAN IN-THE-EAR HEARING AIDS

WAX TRAPS

If the hearing aid sounds weak or there is no amplification, the wax trap may be plugged. To change the wax trap:

1. Remove the white wax trap by placing the empty end of the wax trap tool into the wax trap. Twist slightly and then pull the used wax trap out.
2. Replace the wax trap by flipping the wax trap tool over to the end with the new wax trap. Gently insert the new wax trap into hearing aid, twist slightly, and remove the tool.



THE
SPEECH & HEARING
CENTER

Celebrating 70 Years - 1953-2023

423.622.6900 | www.SpeechHearing.com